



LSO Silicone Oil™



1. Ultra pure ophthalmic grade silicone oil
 - Fractionated to remove low molecular weight silicones
 - No potential for inflammation
 - No emulsification requiring premature removal from the eye
2. Available in 1000cP & 5000cP
3. Available in 10mL vial or 10mL syringe

LSO QuickSet™



- Connects to most VFI machines*
- Fast, easy, convenient & clean dispensing
- No pouring required
- No wastage of oil
- No wasted time waiting for air bubbles to coalesce
- See through holder protects against accidental breakage of glass syringe
- Patented design

LSO QuickSet™ Combo Pack

- LSO Silicone Oil™ and the LSO QuickSet™ can be purchased separately or in the convenient economical LSO QuickSet™ Combo Pack

LSO Silicone Oil™
+
LSO QuickSet™

* Either directly or with the use of an adapter

Labtician Ophthalmics, Inc.
2140 Winston Park Drive,
Unit 6, Oakville, ON,
Canada L6H 5V5
Tel: (905) 829-0055
Fax: (905) 829-0056
1 800 265 8391
e-mail: info@labtician.com
website: www.labtician.com



ISO 13485



ISO 13485 FM72835