

**The healing begins
the moment you begin**



The Perfect Synergy

Monoderm's™ monofilament construction combined with the sharpness of the DermaGlide® needle creates the perfect synergy for smooth passage through the skin. The DermaGlide® tip was designed to minimize trauma, suture line bleeding and scarring.

DermaGlide®



- ▶ Innovative quad-sided diamond geometry
- ▶ Penetrates highly elastic dermal tissue without tenting or drag
- ▶ Holds its fine edge pass after pass

MONODERM™

Monofilament Synthetic
Absorbable Suture

Predictable, Reliable Performance

Monoderm™, Sharpoint's® monofilament synthetic absorbable suture, offers excellent handling properties. Monoderm™ provides minimal resistance during passage through tissue and excellent tensile strength that holds during the critical wound healing period.

- ▶ Pliable, easy-to-tie
- ▶ Elicits lower tissue reaction than gut
- ▶ Predictable strength retention: 64 – 76% at 7 days; 40% at 14 days
- ▶ Fast absorbing – essentially complete at 90 days
- ▶ Frequently used for blepharoplasty, rhinoplasty, tumor removals and laceration repairs

To learn more about the broad line of Sharpoint® suture products, contact our Customer Service Team at:

800 523-3332 or 610 404-1000.

Visit us on the Web at www.sharpoint.com.






Sharpoint®
The Surgeon's Edge



MONODERM™

A monofilament synthetic absorbable suture

DermaGlide® 

NEEDLE	MATERIAL	LENGTH	SIZE			
			6-0	5-0	4-0	3-0
DGL11  11.00mm 3/8 Circle (P-1)	Undyed Monofilament	18"/45cm	YD-489N Y489G	YD-490N Y490G		
DGL13  13.00mm 3/8 Circle (P-3, PC-1)	Undyed Monofilament	18"/45cm	YD-492N Y492G	YD-493N Y493G	YD-494N Y494G	
	Violet Monofilament	18"/45cm		YD-463N Y463G	YD-464N Y464G	
DGL19  19.00mm 3/8 Circle (PS-2, PC-5)	Undyed Monofilament	18"/45cm		YD-495N Y495G	YD-496N Y496G	YD-497N Y497G
	Violet Monofilament	18"/45cm			YD-513N Y513G	
DGL24  24.00mm 3/8 Circle (PS-1)	Undyed Monofilament	18"/45cm			YD-682N Y682H	YD-936N Y936H
HGL9  9.00mm 1/2 Circle (P-2)	Undyed Monofilament	18"/45cm		YD-503N ■	YD-504N ■	

Product numbers are color coded respectively as Sharpoint® (blue) and Ethicon (black).
 N denotes a box of 12.

For more information, contact our Customer Service Team at 800 523-3332.