


UltraPlug[®]

Punctal Plugs



Punctal Plugs

UltraPlug®

Punctal Plugs



**EXTENDED WEAR
SYNTHETIC ABSORBABLE
(PCL)**

Evaluate within 2 months



COLLAGEN

Evaluate within 1 week



SILICONE

Permanent,
easily removable

Extend Your Treatment Options

UltraPlug UltraPlug UltraPlug UltraPlug UltraPlug UltraPlug UltraPlug



The UltraPlug® Dry Eye Management System provides an eye-opening answer for increasing patient awareness of Dry Eye Syndrome:

- Educate your patients
- Test your patients
- Perform the procedure
- Optimize procedure revenue

ULTRAPLUG® SILICONE PUNCTAL PLUGS

Inserted into the punctal opening to block tear drainage through the canaliculus, silicone plugs provide long-term relief of dry eye symptoms, and enable the eye to retain more of its own natural tears:

- Low-profile cap design virtually eliminates foreign body sensation
- Flexibility of cap design facilitates easy removal with jeweler's forceps, if necessary
- Individually packaged and preloaded on a sterile, disposable inserter / dilator
- Manufactured from soft silicone utilizing a single piece mold design for greater patient comfort




PPLUG4	0.4mm plug / inserter (2 / box)
PPLUG5	0.5mm plug / inserter (2 / box)
PPLUG6	0.6mm plug / inserter (2 / box)
PPLUG7	0.7mm plug / inserter (2 / box)
PPLUG8	0.8mm plug / inserter (2 / box)

ULTRAPLUG® EXTENDED WEAR SYNTHETIC ABSORBABLE PUNCTAL PLUGS

Inserted into the punctal opening to block tear drainage through the canaliculus, synthetic absorbable plugs provide extended temporary relief of dry eye symptoms:

- Indicated for use in the treatment of Dry-Eye Syndrome and the dry eye components of various ocular surface diseases
- May be used after ocular surgery to prevent dry eye related complications and to enhance the retention of ocular medications
- Ideal for post-op dry eye management of LASIK patients
- Treatment option for seasonal allergy patients experiencing dry eye symptoms
- May be beneficial for contact lens patients experiencing dry eye related symptoms
- Manufactured from E-Caprolactone-L-Lactide copolymer (PCL)
- Essentially absorbed within 2-6 months



Standard Length - 2.0mm

EPLUG2	0.2mm diameter / 2.0mm length (2 plugs per package / 10 packages per box)
EPLUG3	0.3mm diameter / 2.0mm length (2 plugs per package / 10 packages per box)
EPLUG4	0.4mm diameter / 2.0mm length (2 plugs per package / 10 packages per box)

Short Length - 1.6mm


ESPLUG2	0.2mm diameter / 1.6mm length (2 plugs per package / 10 packages per box)
ESPLUG3	0.3mm diameter / 1.6mm length (2 plugs per package / 10 packages per box)
ESPLUG4	0.4mm diameter / 1.6mm length (2 plugs per package / 10 packages per box)

If you would like to place an order or require assistance cross-referencing one of our products, please call 800 523-3332 or 610 404-1000.

ULTRAPLUG® COLLAGEN ABSORBABLE PUNCTAL PLUGS

Inserted into the punctal opening to block tear drainage through the canaliculus, collagen plugs provide temporary relief of dry eye symptoms:

- Manufactured from precisely measured premium collagen
- Used as a diagnostic assessment for permanent punctal occlusion
- Each individually packaged card contains 6 plugs




Standard Length - 2.0mm

C-PLUG2	0.2mm diameter / 2.0mm length (6 plugs per package / 10 packages per box)
C-PLUG3	0.3mm diameter / 2.0mm length (6 plugs per package / 10 packages per box)
C-PLUG4	0.4mm diameter / 2.0mm length (6 plugs per package / 10 packages per box)
C-PLUG-V	4 cards each of 3 different sizes - 0.2, 0.3, 0.4mm diameter / 2.0mm length (6 plugs per package / 12 packages per box)

ULTRAPLUG® EVALUATION / INSERTION PRODUCTS

- A stepped system of Gauge / Dilator instruments provides a convenient and controlled method for punctal measurement and dilation. Available as a non-sterile autoclavable kit or as non-sterile individual instruments
- Recommended for use exclusively with UltraPlug® Punctal Plugs



PAK-101

UltraPlug® Punctal Gauge / Dilator Set: a set of 3 instruments, 2 sizes on each instrument. First instrument: 0.4mm / 0.5mm, second instrument: 0.6mm / 0.7mm, third instrument: 0.8mm / 0.9mm. Packaged non-sterile in an autoclavable case



PGAUGE4/5

0.4mm / 0.5mm UltraPlug® Punctal Gauge / Dilator

PGAUGE6/7

0.6mm / 0.7mm UltraPlug® Punctal Gauge / Dilator

PGAUGE8/9

0.8mm / 0.9mm UltraPlug® Punctal Gauge / Dilator

- Zone-Quick is an easy and accurate system that measures tear volume in just 15 seconds per eye. One kit can test up to 50 patients
- Easy insertion and removal
- Packaged 50 sterile sets to a box / 2 threads per sets



Z-QUICK

Phenol red thread tear test strips



ULTRAPLUG® PUNCTAL PLUG DRY EYE MANAGEMENT SYSTEM (Not available in all countries)

System Contains:

Silicone Plugs:

- 3 Boxes - 0.4mm Plug
- 4 Boxes - 0.5mm Plug
- 2 Boxes - 0.6mm Plug
- 1 Box - 0.7mm Plug

Collagen Plugs:

- 1 Box - Variety Pack, Standard Length (0.2, 0.3, 0.4mm diameters)

Extended Wear Synthetic Absorbable Plugs:

- 1 Box - Standard Length (0.3mm diameter)

Gauging & Testing Materials:

- 1 Box - Zone-Quick (tear test strips)
- 1 Gauge / Dilator Set

Educational / Promotional Materials:

- 1 UltraPlug® Punctal Plug Video
- 1 UltraPlug® Punctal Plug Poster
- 1 UltraPlug® Punctal Plug Acrylic Display
- 50 UltraPlug® Punctal Plug Patient Pamphlets
- 1 UltraPlug® Punctal Plug Medicare Reimbursement / Promotional Information CD

DEKIT101

DEKIT201

<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

If you would like to place an order or require assistance cross-referencing one of our products, please call 800 523-3332 or 610 404-1000.

ULTRAPLUG® PUNCTAL PLUG EDUCATIONAL MATERIALS

UltraPlug® educational materials for staff and patients in a variety of media



OPG-PM-045R1

PHYSICIAN VIDEO - UltraPlug® Punctal Plug - VHS

OPG-PM-050

PHYSICIAN VIDEO - UltraPlug® Punctal Plug - DVD

OPG-PM-172

PHYSICIAN VIDEO - UltraPlug® Punctal Plug - MPEG



OPG-PM-038

UltraPlug® Punctal Plug Patient Pamphlet (50 per pack)

OPG-PM-092


UltraPlug® Punctal Plug Acrylic Pamphlet Holder.
Freestanding display supplied with 8 1/2" x 11" color insert.
Patient Pamphlets ordered separately



OPG-PM-091

UltraPlug® Punctal Plug Patient Poster (16" x 20")

NOTES:



UltraPlug *UltraPlug* *UltraPlug* *UltraPlug* *UltraPlug* *UltraPlug* *Ultra*

UltraPlug[®] ***Punctal Plugs***



Surgical Specialties Corporation, Ophthalmic Products Group, 100 Dennis Drive, Reading, PA 19606, U.S.A.
610 404-1000 • 800 523-3332 (U.S.A. only) • Fax: 610 404-4010 • www.angiotech.com • OPG-PM-147R2 04/07

UltraPlug® and the Lasik Plug™ are trademarks of Surgical Specialties Corporation.

©2007 Angiotech Pharmaceuticals, Inc. ©2007 Surgical Specialties Corporation. All Rights Reserved.

UltraPlug®
Punctal Plugs