



**QUILL SRS: ORDER FORM**

<b>Surgery name</b>			<b>Date</b>		
<b>Address</b>			<b>Contact Name</b>		
<b>Phone Number</b>			<b>Order #</b>		
<b>Fax Number</b>					

QUILL CODE	SIZE	NEEDLE SIZE	SUTURE LENGTH	BOX SIZE	QTY
------------	------	-------------	---------------	----------	-----

**ABSORPBABLE QUILL SRS**

**MONODERM, Diamond Point      ABSORBABLE UNDYED      90-120 days total absorption**  
**Tensile Strength 62% 7 days post implantation, 27% Tensile Strength 14 days post implantation**

YA-1019Q	3-0	18mm 3/8 Circle	3.5cm x 3.5cm	12	
YA-1018Q	2-0	18mm 3/8 Circle	3.5cm x 3.5cm	12	
YA-1001Q	3-0	18mm 3/8 Circle	7cm x 7cm	12	
YA-1000Q	2-0	18mm 3/8 Circle	7cm x 7cm	12	
YA-1016Q	3-0	18mm 3/8 Circle	14cm x 14cm	12	
YA-1012Q	2-0	18mm 3/8 Circle	14cm x 14cm	12	
YA-1014Q	0	18mm 3/8 Circle	14cm x 14cm	12	
YA-1011Q	2-0	26mm 3/8 Circle	7cm x 7cm	12	
YA-1010Q	0	26mm 3/8 Circle	7cm x 7cm	12	
YA-1015Q	2-0	26mm 3/8 Circle	14cm x 14cm	12	
YA-1014Q	0	26mm 1/2 Circle	14cm x 14cm	12	

**MONODERM, Taper Point      ABSORBABLE UNDYED      90-120 days total absorption**  
**Tensile Strength 62% 7 days post implantation, 27% Tensile Strength 14 days post implantation**

YA-1003Q	0	36mm 1/2 Circle	14cm x 14cm	12	
----------	---	-----------------	-------------	----	--

**PDO (Polydioxanone), Diamond Point      ABSORBABLE DYED (Violet #2)      180 days total absorption**  
**Tensile Strength post implantation 3-0 & larger - 80% @ 2 weeks; 80% @ 4 weeks; 40-70% @ 6 weeks**

RA-1017Q	3-0	12mm 3/8 Circle	7cm x 7cm	12	
RA-1008Q	3-0	18mm 3/8 Circle	3.5cm x 3.5cm	12	
RA-1002Q	2-0	18mm 3/8 Circle	3.5cm x 3.5cm	12	
RA-1010Q	3-0	18mm 3/8 Circle	7cm x 7cm	12	
RA-1007Q	2-0	18mm 3/8 Circle	7cm x 7cm	12	
RA-1050Q	3-0	18mm 3/8 Circle	14cm x 14cm	12	
RA-1046Q	2-0	18mm 3/8 Circle	14 cm x 14cm	12	
RA-1051Q	3-0	18mm 3/8 Circle	24cm x 24cm	12	
RA-1036Q	2-0	18mm 3/8 Circle	24cm x 24cm	12	
RA-1016Q	2-0	26mm 3/8 Circle	7cm x 7cm	12	
RA-1015Q	0	26mm 3/8 Circle	7cm x 7cm	12	
RA-1021Q	3-0	26mm 3/8 Circle	14cm x 14cm	12	
RA-1020Q	2-0	26mm 3/8 Circle	14cm x 14cm	12	
RA-1019Q	0	26mm 3/8 Circle	14cm x 14cm	12	
RA-1052Q	3-0	26mm 3/8 Circle	24cm x 24cm	12	
RA-1028Q	2-0	26mm 3/8 Circle	24cm x 24cm	12	
RA-1030Q	0	26mm 3/8 Circle	24cm x 24cm	12	
RA-1001Q	0	18mm 1/2 Circle	14cm x 14cm	12	
RA-1004Q	1	26mm 1/2 Circle	14cm x 14cm	12	
RA-1006Q	2	26mm 1/2 Circle	14cm x 14cm	12	
RA-1056Q	3-0	26mm 1/2 Circle	24cm x 24cm	12	
RA-1032Q	1	26mm 1/2 Circle	24cm x 24cm	12	
RA-1034Q	2	26mm 1/2 Circle	24cm x 24cm	12	

<b>PDO (Polydioxanone), Taper Cutting    ABSORBABLE DYED (Violet #2)    180 days total absorption</b>					
<b>Tensile Strength post implantation 3-0 &amp; larger - 80% @ 2 weeks; 80% @ 4 weeks; 40-70% @ 6 weeks</b>					
RA-1013Q	2-0	50mm Straight	7cm x 7cm	12	
RA-1023Q	0	50mm Straight	7cm x 7cm	12	
RA-1024Q	1	50mm Straight	7cm x 7cm	12	
<b>PDO (Polydioxanone), Taper Point    ABSORBABLE DYED (Violet #2)    180 days total absorption</b>					
<b>Tensile Strength post implantation 3-0 &amp; larger - 80% @ 2 weeks; 80% @ 4 weeks; 40-70% @ 6 weeks</b>					
RA-1000Q	2-0	36mm 1/2 Circle	14cm x 14cm	12	
RA-1003Q	1	36mm 1/2 Circle	14cm x 14cm	12	
RA-1005Q	2	36mm 1/2 Circle	14cm x 14cm	12	
RA-1029Q	0	36mm 1/2 Circle	24cm x 24cm	12	
RA-1031Q	1	36mm 1/2 Circle	24cm x 24cm	12	
RA-1033Q	2	36mm 1/2 Circle	24cm x 24cm	12	
<b>NON ABSORBABLE QUILL SRS</b>					
<b>NYLON, Diamond Point    NON-ABSORPBABLE (FD&amp;C Blue #2)</b>					
AA-1003Q	1	26mm 1/2 Circle	14cm x 14cm	12	
AA-1005Q	2	26mm 1/2 Circle	14cm x 14cm	12	
<b>POLYPROPYLENE, Diamond Point    NON-ABSORPBABLE (Clear)</b>					
JA-1006Q	2-0	18mm 3/8 Circle	7cm x 7cm	12	
JA-1005Q	1	26mm 1/2 Circle	14cm x 14cm	12	
JA-1001Q	2	26mm 1/2 Circle	14cm x 14cm	12	

**PRICES EXCLUDES GST**