

WOUND CLOSURE CATALOGUE

Full range of suture material for all of your wound closure needs



High-quality, consistently sharp needles.

- ▶ Extra sharp cutting edges extending from the point to its widest portion reduce tissue drag.
- ▶ Precise, smooth penetration through tissue minimizes trauma.
- ▶ Needle channel is fitted to the specific size of the suture material.

DermaGlide™ Needles

Innovative quad sided diamond geometry provides unsurpassed incision and healing characteristics for enhanced cosmesis.

| SharpPoint | Ethicon | USS | |
|-------------------|-----------|------|------|
| 3/8 CIRCLE | | | |
| DGL11 | P-1 | P-10 | 11mm |
| DGL13 | P-3/PC-1 | ■ | 13mm |
| DGL19 | PS-2/PC-5 | P-12 | 19mm |

Precision Reverse Cutting Needles

Precision sharpened premium quality needle.

| SharpPoint | LOOK | Ethicon | USS | |
|-------------------|------|---------|------|------|
| 3/8 CIRCLE | | | | |
| DSM11 | C1 | P-1 | P-10 | 11mm |
| DSM13 | C3 | P-3 | ■ | 13mm |
| DSM16 | PC34 | PS-3 | P-11 | 16mm |
| DSM18 | PC31 | PS-2 | P-12 | 18mm |
| DSM24 | PC33 | PS-1 | P-14 | 24mm |

UltraGlide® Needles

UltraGlide™ is a SharpPoint™ patented needle design. Utilises our Infinite Edge™ technology sharpening process. Modified spatula needle facilitates easier knot burial for corneal grafts.

| SharpPoint | Ethicon | USS | |
|-------------------|---------|----------|--------|
| 3/8 CIRCLE | | | |
| DGL6 | CS140-6 | SE-140-6 | 6.15mm |
| DGL7 | ■ | SE-140-6 | 6.60mm |

Reverse Cutting Needles

| SharpPoint | LOOK | Ethicon | USS | |
|-------------------|------|---------|------|------|
| 3/8 CIRCLE | | | | |
| DS14 | ■ | ■ | ■ | 14mm |
| DS16 | C22 | FS-3 | C-12 | 16mm |
| DS18/Fine | C6 | FS-2 | C-13 | 18mm |
| DS24 | C7 | FS-1 | C-14 | 24mm |
| C8/DS26 | C8 | FS | C-15 | 26mm |



ABSORBABLE

NON ABSORBABLE

Selected products from a comprehensive range of sutures

SHARPOINT LOOK ETHICON SIZE NEEDLE SUTURE SINGLE BOX
FSSB CODE CODE SIZE LENGTH / DOUBLE ARMED SIZE

| PCL ABSORBABLE SUTURE MATERIAL | | | | | | |
|---|-------------|-----|------|------|----|----|
| MONOFILAMENT SYNTHETIC - ABSORBABLE UNDYED | | | | | | |
| Extended Wound Support 77% of Strength at 21 Days | | | | | | |
| DERMAGLIDE NEEDLE 3/8 CIRCLE | | | | | | |
| QD-489N | Z489G | 6-0 | 11MM | 45CM | SA | 12 |
| QD-493N | Z493 | 5-0 | 13MM | 45CM | SA | 12 |
| QD-494N | Z494G | 4-0 | 13MM | 45CM | SA | 12 |
| QD-495N | | 5-0 | 19MM | 45CM | SA | 12 |
| QD-496N | Z496G/Z823G | 4-0 | 19MM | 45CM | SA | 12 |
| QD-497N | Z497G/Z824G | 3-0 | 19MM | 45CM | SA | 12 |
| QD-682N | Z682G | 4-0 | 24MM | 45CM | SA | 12 |

MONODERM MONOFILAMENT SYNTHETIC ABSORBABLE UNDYED
High Initial Tensile Strength Diminishing to 20-30% at 14 Days

| DERMAGLIDE NEEDLE 3/8 CIRCLE | | | | | | |
|------------------------------|-------|-----|------|------|----|----|
| YD-489N | Y489G | 6-0 | 11MM | 45CM | SA | 12 |
| YD-490N | Y490G | 5-0 | 11MM | 45CM | SA | 12 |
| YD-463N Violet | Y463G | 5-0 | 13MM | 45CM | SA | 12 |
| YD-492N | Y492G | 6-0 | 13MM | 45CM | SA | 12 |
| YD-493N | Y493G | 5-0 | 13MM | 45CM | SA | 12 |
| YD-494N | Y494G | 4-0 | 13MM | 45CM | SA | 12 |
| YD-495N | Y495G | 5-0 | 19MM | 45CM | SA | 12 |
| YD-496N | Y496G | 4-0 | 19MM | 45CM | SA | 12 |
| YD-497N | Y497G | 3-0 | 19MM | 45CM | SA | 12 |

| PRECISION REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|---|--|-----|------|------|----|----|
| YP-489N | | 6-0 | 11MM | 45CM | SA | 12 |
| YP-463N | | 5-0 | 13MM | 45CM | SA | 12 |
| YP-493N | | 5-0 | 13MM | 45CM | SA | 12 |
| YP-492N | | 6-0 | 13MM | 45CM | SA | 12 |

POLYGLYCOLIC ACID (PGA) - BRAIDED, SYNTHETIC ABSORBABLE
(Similar to Vicryl) 50% of Strength at 14 Days

| FSSB SPATULA NEEDLE 1/4 CIRCLE | | | | | | |
|--------------------------------|-------|-----|-------|------|----|----|
| A007PGA | J571G | 5-0 | 8.5MM | 45CM | DA | 12 |
| A001PGA | J570G | 6-0 | 8.0MM | 45CM | DA | 12 |

| FSSB SPATULA NEEDLE 1/2 CIRCLE | | | | | | |
|--------------------------------|--|-----|-------|------|----|----|
| A069PGA | | 6-0 | 6.4MM | 30CM | DA | 12 |

| FSSB SPATULA NEEDLE 3/8 CIRCLE | | | | | | |
|--------------------------------|-------|------|-------|------|----|----|
| A034PGA | J576 | 7-0 | 6.4MM | 45CM | DA | 12 |
| A003PGA | J547G | 8-0 | 6.4MM | 30CM | DA | 12 |
| A063PGA | V549G | 9-0 | 6.4MM | 30CM | SA | 12 |
| A002PGA | V960G | 10-0 | 6.2MM | 10CM | SA | 12 |

| FSSB TAPER POINT NEEDLE 3/8 CIRCLE | | | | | | |
|------------------------------------|-------|-----|-------|------|----|----|
| E007PGA | J401G | 8-0 | 6.4MM | 30CM | SA | 12 |

| PRECISION REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|---|-------|-----|------|------|----|----|
| LOOK-492B | J492G | 6-0 | 13MM | 45CM | SA | 12 |
| LOOK-493B | J493G | 5-0 | 13MM | 45CM | SA | 12 |
| LOOK-494B | J494G | 4-0 | 13MM | 45CM | SA | 12 |
| LOOK-496B | J496G | 4-0 | 19MM | 45CM | SA | 12 |
| LOOK-497B | J497G | 3-0 | 19MM | 45CM | SA | 12 |

| REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|-----------------------------------|-------|-----|------|------|----|----|
| LOOK-421B | J421H | 5-0 | 19MM | 45CM | SA | 12 |
| LOOK-422B | J422H | 4-0 | 19MM | 45CM | SA | 12 |
| LOOK-423B | J423H | 3-0 | 24MM | 45CM | SA | 12 |
| LOOK-441B | J441H | 4-0 | 24MM | 45CM | SA | 12 |
| LOOK-442B | J423H | 3-0 | 24MM | 70CM | SA | 12 |
| LOOK-443B | J422H | 2-0 | 24MM | 70CM | SA | 12 |

| SKIN ADHESIVE | | | |
|----------------------|--|--|------|
| SS 035T | SURGISEAL TOPICAL SKIN ADHESIVE 0.35ML | | EACH |
| SS 035TX | SURGISEAL TOPICAL SKIN ADHESIVE 0.35ML | | 12 |

| SKIN MARKERS | | | |
|---------------------|-------------------------|--|----|
| 52-100 | GENTIAN INK MARKER PENS | | 10 |

SHARPOINT LOOK ETHICON SIZE NEEDLE SUTURE SINGLE BOX
FSSB CODE CODE SIZE LENGTH / DOUBLE ARMED SIZE

| NYLON NON ABSORBABLE SUTURE MATERIAL | | | | | | |
|---|--|------|-----|------|----|----|
| BLACK MONOFILAMENT | | | | | | |
| M.E.T POINT NEEDLE 3/8 CIRCLE | | | | | | |
| AA-0111 | | 10-0 | 6MM | 13CM | SA | 12 |
| AA-0121 | | 9-0 | 6MM | 13CM | SA | 12 |
| AA-0122 | | 8-0 | 6MM | 13CM | SA | 12 |

| DERMAGLIDE NEEDLE 3/8 CIRCLE | | | | | | |
|------------------------------|-------------|-----|------|------|----|----|
| AD-1698 | 1677G/1894G | 6-0 | 13MM | 45CM | SA | 12 |
| AD-1698S | | 6-0 | 13MM | 20CM | SA | 12 |
| AD-1666 | 1666G/1895G | 5-0 | 19MM | 45CM | SA | 12 |

| ULTRAGLIDE NEEDLE 3/8 CIRCLE | | | | | | |
|------------------------------|-------------|------|-----|------|----|----|
| AA-2627N | 9003G W1770 | 10-0 | 6MM | 30CM | DA | 12 |
| AA-2673N | | 10-0 | 6MM | 10CM | SA | 12 |

| ULTRAGLIDE NEEDLE 1/2 CIRCLE | | | | | | |
|------------------------------|-------|------|-----|------|----|----|
| AA-2713N | 9016G | 11-0 | 5MM | 15CM | DA | 12 |

| PRECISION REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|---|------|-----|------|------|--|--|
| AA-2177N | 698G | 5-0 | 13MM | 45CM | | |

| REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|-----------------------------------|------|-----|------|------|----|----|
| AC-0525D | 667G | 6-0 | 12MM | 45CM | SA | 36 |
| AC-0525N | 667G | 6-0 | 18MM | 45CM | SA | 12 |
| AC-0568D | 661G | 5-0 | 18MM | 45CM | SA | 36 |
| AC-0568N | 661G | 5-0 | 18MM | 45CM | SA | 12 |
| AC-0182N | 622G | 4-0 | 18MM | 45CM | SA | 12 |
| AC-0182D | 662G | 4-0 | 18MM | 45CM | SA | 36 |
| AC-0586N | 669G | 3-0 | 24MM | 75CM | SA | 12 |
| AC-0586D | 669G | 3-0 | 24MM | 75CM | SA | 36 |

POLYPROPYLENE BLUE MONOFILAMENT
Inert, Smooth, Better Knot Security

| DERMAGLIDE NEEDLE 3/8 CIRCLE | | | | | | |
|------------------------------|-------------|-----|------|------|----|----|
| JD-8697 | 8997G | 6-0 | 11MM | 45CM | SA | 12 |
| JD-8694N | | 5-0 | 11MM | 45CM | SA | 12 |
| JD-8698 | 8698G/8618G | 5-0 | 13MM | 45CM | SA | 12 |
| JD-8686 | 8686G/8630G | 5-0 | 19MM | 45CM | SA | 12 |

| PRECISION REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|---|-------|-----|------|------|----|----|
| LOOK-1082B | 8695G | 6-0 | 13MM | 45CM | SA | 12 |

| REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|-----------------------------------|-------|-----|------|------|----|----|
| LOOK-1011B | | 6-0 | 12MM | 25CM | SA | 12 |
| LOOK-1018B | 8661G | 5-0 | 18MM | 45CM | SA | 12 |
| JC-1227D | 8661G | 5-0 | 18MM | 45CM | SA | 36 |
| JC-1229N | 8683G | 4-0 | 18MM | 45CM | SA | 12 |
| JC-1229D | 8683G | 4-0 | 18MM | 45CM | SA | 36 |
| JC-1231N | | 3-0 | 18MM | 45CM | SA | 12 |
| JC-1231D | | 3-0 | 18MM | 45CM | SA | 36 |
| JC-1235N | 8684G | 3-0 | 24MM | 45CM | SA | 12 |
| JC-1239N | | 3-0 | 26MM | 45CM | SA | 12 |

SILK BLACK BRAIDED
Excellent Handling

| PRECISION REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|---|------|-----|------|------|----|----|
| DA-2622N | | 6-0 | 7MM | 45CM | DA | 12 |
| DA-0187N | 639G | 6-0 | 11MM | 45CM | SA | 12 |
| DA-0188N | | 6-0 | 11MM | 45CM | DA | 12 |
| LOOK-1265B | | 5-0 | 13MM | 45CM | SA | 12 |

| REVERSE CUTTING NEEDLE 3/8 CIRCLE | | | | | | |
|-----------------------------------|-------|-----|------|------|----|----|
| DC-0210N | 683GG | 4-0 | 18MM | 45CM | SA | 12 |
| DC-0156N | 684G | 3-0 | 24MM | 45CM | SA | 12 |
| DC-0212N | | 3-0 | 26MM | 75CM | SA | 12 |
| DC-0211N | 685G | 2-0 | 26MM | 45CM | SA | 12 |

POLYESTER/POLYVIOLENE BRAIDED
High Strength for Long Term Support

| SPATULA NEEDLE 1/2 CIRCLE | | | | | | |
|---------------------------|-------|-----|-----|------|----|----|
| A003PEW White | 1760G | 5-0 | 8MM | 45CM | DA | 12 |

N = BOX 12 SUTURES
D = BOX 36 SUTURES

Critical Assist are exclusive suppliers of **Quill** Self Retaining Systems.

Contact one of our team to discuss your clinical requirements