


## QRS UNIVERSAL ECG WITH OFFICE MEDIC PATIENT SOFTWARE



- + Highly portable & affordable PC based ECG.
- + 12 lead interpretive universal ECG allows practitioners to perform important diagnostic tests from a lap top, PC or pocket PC.
- + Automatic narrative interpretation and measurement analysis using the advanced Louvain Algorithm, which has the **best total accuracy** when compared to leading competitors.
- + Print full page reports on A4 paper.
- + Review and zoom into data with electronic callipers, enter comments and modify the interpretation before confirming the results.
- + Optional Bluetooth wireless device available.



## RADIAL PULSE WAVE ACQUISITION DEVICE BPRO® 24HR Ambulatory Blood Pressure Monitoring

- + Does not disrupt sleep
- + Wearer is not aware that BP readings are being taken
- + No inflatable cuff to be worn



### A more complete assessment of the hypertensive condition

- + 24-hour ambulatory BP patterns (macroscopic analysis)
- + Arterial pulse waveform analysis (microscopic analysis)

### An evidence-based approach to diagnosis & management of hypertension

- + Identifies BP patterns.
- + Enables individualized treatment

More information on  
[www.healthstats.com](http://www.healthstats.com)

## FUKUDA DENSHI CARDIMAX ECG's Small light-weight & portable ECG Monitors

**FX 7102** Non Interpretative Model      **FCP 7101** Interpretive Model

- + Large LCD display offers unparalleled ease of use for checking patient data.
- + Clear, accurate viewing of ECG waveforms.
- + Arrhythmia analysis and interpretation is a standard feature on the FCP-7101.
- + Extended battery life up to 3 Hours internal Memory for 128 ECG
- + Built-in large LCD (320 x 240 dots)
- + 1/3 channel on 63mm paper.
- + LAN support for external communication

CONSUMABLES AVAILABLE



FX 7102



FCP 7101

### FX 7402

- + Time-saving A4 size recording format
- + Clear ECG display
- + Touch Screen for easy operation
- + Alpha-numeric keyboard
- + Interpretation program option
- + Two way power supply AC and rechargeable battery (option).
- + Stress Test (option).
- + Stress/Interpretative cartridges additional cost



### Electrodes 4530M

- + Unirest 'All in One' electrodes are pre-gelled.
- + Connect directly to the ECG lead and eliminate the need for alligator clips.
- + Less patient prep required.

## EXERGEN TEMPORAL ARTERY THERMOMETRY

### TAT2000 & TAT 5000A (Hospital model)

- + Latest infrared technology seeks out & locks onto the highest temperature when swept along temporal artery and touched behind ear lobe.
- + A built in thermistor measures the ambient temperature and computer synthesizes both to find the core temperature. Ideal for use in ED.
- + Extremely accurate, reliable and robust. Minimal downtime.
- + More accurate than ear thermometry (Harvard Medical School study).
- + Non invasive, especially good for babies, children and elderly patients.
- + Silver Antimicrobial head. Common microbes killed instantly, VRE & MRSA undetectable in 6 hours.
- + Huge cost saving attached with probe covers being unnecessary for majority of patients. Just wipe with alcohol wipes between patients.
- + Real time measurement like a pulmonary artery thermometry and as accurate.



**COMPREHENSIVE RANGE OF SUTURES AND SURGICAL INSTRUMENTS ALSO AVAILABLE**

## HADECO MEDICAL ULTRASONIC DOPPLERS

### SMAES101EX2 Echo Sounder Hand-held Fetal Doppler

- + Basic cell-phone sized fetal Doppler.
- + Built-in speaker and transducer.
- + Detects fetal heart beat at 10-12 weeks (non-waterproof).
- + Carry case and gel included
- + One-handed operation.
- + Large LCD for numerical data, PEAK & MEAN velocities.
- + Freeze key freezes heart rate for notation.
- + Safe and effective low-intensity beam.
- + Fixed 2Mhz Fetal probe



- + Heart rate 60-240 BPM
- + Uses 2 AA batteries

### SMAES100V11 Obstetric Minidop

- + 2MHz obstetric Doppler
- + water-proof probe.
- + Fetal pocket Doppler with 80% improved great sounds.
- + Displays fetal heart rate with LCD & speaker.
- + 3 modes available: Average, Manual counting, & Real-time.
- + Safe and effective low intensity beam,
- + Freeze key for heart rate notation
- + 9 volt Low battery indicator



### SMAES100V3 Bi-directional Arterial Doppler

- + Large backlit LCD screen displays real-time waveform & numerical data
- + Advanced easy operating control. Headphone for operation in noisy environments.
- + 30 waveform memory and USB interface.
- + Connects to a PC with additional software (additional price).
- + Uses a 9 volt battery
- + Can be used as a Fetal Doppler with a 2MHz probe (not water-proof).
- + Carry case & gel included.



### SMASD45 Arterial Minidop

- + Top- of the line hand-held, portable arterial Doppler.
- + Integrated printer.
- + Superior sound quality with 40% improved great sounds.
- + Large LCD display showing wave form and numerical parameters.
- + Built in lithium re-chargeable battery.
- + 30 waveform memory and USB interface.
- + PPG21 infrared probe available (additional price).



## ISHIHARA COLOUR VISION TEST

Two editions available:

- + ISHI 14 Plate for quick screening.
- + ISHI 24 plate for detailed investigation



## QRS SPIROMETRY – The only PC Card Spirometer in the world.

### Z-7000-0100 SpiroCard®

The SpiroCard brings full function Spirometry to computers and handhelds.

- + Smallest, lightest and most accurate Spirometer on the market.
- + Perform full-function Spirometry anywhere!
- + Connect via PC Card slot or USB



### Z-7000-0200 SpirOxCard®

The first combination Spirometer / Pulse Oximeter designed for the mobile healthcare professional.

- + Perform diagnostic testing for asthma COPD and other respiratory diseases in seconds.
- + Unparalleled accuracy.
- + Connect via PC Card Slot or USB.



### QRS Pre-Calibrated Mouthpieces Z-5000-2608

Patented pre-calibrated mouthpieces eliminate the need for sterilization and ongoing calibration.

## CONSUMABLES AVAILABLE

## PULSE OXIMETERS

### Critical Assist Handheld Pulse Oximeter – SA210

For spot checking or attended-care monitoring of adult or pediatric patients.

- + Ideal for sleep apnoea screening – non invasive
- + Motion artifact
- + Lightweight and easy to use
- + Rolling 288 data memory
- + Oxygen Sat. Range: 35%-100%
- + Accuracy: 70-100%, ± 2 digital; 35-69%, ± 3 digital
- + Pulse Rate Range: 30-250 BPM
- + Accuracy: ± 3 digital
- + Digital LCDs for SPO2 and pulse rate
- + Optional paediatric and neonatal probes.
- + USB port.



### Critical Assist Fingertip Pulse Oximeter – SB100

Health at your fingertips!

- + SpO2 reads oxygen concentration non-invasively.
- + Economical and portable.
- + Heart rate
- + Oxygen concentration
- + Digital LCD's for SPO2 and pulse rate



**COMPREHENSIVE RANGE OF SUTURES AND SURGICAL INSTRUMENTS ALSO AVAILABLE**